

REBOTEC

AUGSBURG 24" SELF PROPELLED Shower commode chair for children



- Custom design for children for use in the shower, over the toilet and at bedside.
- Soft seat and backrest
- Upholstered seat with hygiene opening
- Height adjustable, swinging and detachable foot support
- Swinging armrests for easy side transfer
- 24" tubeless tires

MPN

Augsburg 4" Castors
339.04.97

Augsburg height adjustable head-rest
444.03.39

Augsburg seat width and depth adjustment kit
444.00.97

Augsburg 24" Self Propelled
339.24.97

ACCESSORIES

Bucket Pail
301.00

Safety Belt
322.30.10

IV Pole
303.62

Armrest Cushion
440.10.32

CHAIR DEPTH, WIDTH 97x66 cm

SEAT DEPTH, WIDTH, HEIGHT 35x40x55 cm

TOILET HEIGHT WITHOUT BUCKET DEVICE 48 cm

TOILET HEIGHT WITH BUCKET DEVICE 44 cm

TOILET WIDTH 32 cm

MAX BODY WEIGHT 80 kg



Australian Government
Department of Health
Therapeutic Goods Administration

ARTG 297301

info@rebotec.com.au
rebotec.com.au
218/354 Eastern Valley Way
Chatswood 2067 NSW

