

REBOTEC

CRUTCH PADS

Upholstered Crutch Handle Pads



- Crutch handle upholstered pad
- Pressure-reducing and skin-friendly material
- Slip-resistant neoprene with nylon cover
- Easy Velcro attachment around crutch handle
- Hand wash in warm water
- Sold in pairs

MPN

Small
For Kids and Light crutches
160.01.01

Standard
For crutches with Standard handle
160.01.02

Large
For BIG 150 and BIG 250
160.01.03

Anatom
For Anatom crutches
160.01.04



Australian Government
Department of Health
Therapeutic Goods Administration

ARTG 297302

info@rebotec.com.au
rebotec.com.au
218/354 Eastern Valley Way
Chatswood 2067 NSW

