

REBOTEC

QUADRO

Stable quad forearm crutch with open cuff



- Compression fitting for strong and confident stride without weak weld points
- Wide and stable footprint with 4 interdependent rubber ferrules
- Flexible height with 10 adjustment settings

MPN

Light Grey
155.00

Black
155.10

Accessories

Handle Pads
160.01.02x2



Australian Government
Department of Health
Therapeutic Goods Administration

ARTG 297302

info@rebotec.com.au
rebotec.com.au
218/354 Eastern Valley Way
Chatswood 2067 NSW

HEIGHT FROM HANDLE TO CUFF 24 – 27 cm

HEIGHT OF HANDLE 67 – 89 cm

HEIGHT ADJUSTMENTS x10

