

REBOTEC

QUADRO CANE

Stable quad walking cane



Excellent for users requiring a higher level of support

- Compression fitting for strong and confident stride without weak weld points and minimum rattle
- Wide and stable footprint with 4 interdependent rubber ferrules
- Flexible height with 10 adjustment settings
- 100kg Capacity

MPN

Black
153.10

Light Grey
153.00



Australian Government
Department of Health
Therapeutic Goods Administration

ARTG 297302

info@rebotec.com.au
rebotec.com.au
218/354 Eastern Valley Way
Chatswood 2067 NSW

HEIGHT	63 - 86 cm
--------	------------

DIAMETER Ø IN MM:	19mm
-------------------	------

HEIGHT ADJUSTMENTS	x10
--------------------	-----

