

# REBOTEC

## BIG 250

### Heavy duty forearm crutches



- Heavy duty forearm crutches
- Certified for up to 250kg support with locking system
- Super strong, lightweight aluminium alloy tube
- 10 height adjustments and 3 cuff length adjustments

#### MPN

BIG 250  
106.00.10

BIG 250 Tall  
106.04.10

#### ACCESSORIES

Upholstered Crutch Handle Pad  
160.01.03

Big Foot Ferrules, 22-24mm  
168.10.22

Crutch Lanyard, Charcoal  
167.00.16

	Regular	Tall
HEIGHT FROM HANDLE TO CUFF	22-25cm	22-25cm
HEIGHT OF HANDLE	74-97cm	84-107cm
DIAMETER Ø IN MM	25/22	25/22
MAX BODY WEIGHT	250kg	250kg



Australian Government  
Department of Health  
Therapeutic Goods Administration

ARTG 297302

info@rebotec.com.au  
rebotec.com.au  
218/354 Eastern Valley Way  
Chatswood 2067 NSW

