REBOTEC

BREMEN

Shower commode, PU seat for small adults



- PU cushioned seat, no cracking or reupholstering
- Dignity hygiene aperture cover included
- Commode pail bucket included
- Swing up armrests for easy side transfer

MPN

Blue 357.54.10

ACCESSORIES

Bucket Pail 301.00

Safety Belt 322.30.10

Armrest Cushion 440.10.32

WIDTH, DEPTH	52cm, 95cm
SEAT WIDTH, DEPTH, HEIGHT	42cm, 44cm, 52cm
MAX. TOILET WIDTH	33cm
TOILET HEIGHT WITH BUCKET RAILS	40cm
NET WEIGHT	16.4kg
MAX BODY WEIGHT	90kg



info@rebotec.com.au rebotec.com.au 218/354 Eastern Valley Way Chatswood 2067 NSW

