

# REBOTEC

## BREST

### Over toilet seat frame chair with armrests



- Flexible design is suitable as a toilet seat frame or shower stool.
- Portable Lightweight design folds quickly and easily for travelling
- Height adjustable with sturdy arm-rests for enhanced user security & safety

#### MPN

Brest  
362.00.20

#### ACCESSORIES

Chair Travel Bag  
284.00.30

SEAT HEIGHT 42 cm – 57 cm

SEAT WIDTH 46 cm

SEAT DEPTH 40 cm

REQUIRED SPACE (FOOTPRINT) 45 cm x 59 cm

NET WEIGHT 3.8 kg

MAX BODY WEIGHT 130 kg



Australian Government  
Department of Health  
Therapeutic Goods Administration

ARTG 297301

info@rebotec.com.au  
rebotec.com.au  
218/354 Eastern Valley Way  
Chatswood 2067 NSW

