

REBOTEC

CANNES

Shower stool with armrests



Modular shower stool with armrests and front hygiene opening.

- Stable, assistive support while showering
- Height adjustable with front cut out enabling easy personal hygiene access
- Removable clip in place arms and legs allow for easy travelling

MPN

Cannes
363.00.20

ACCESSORIES

Chair Travel Bag
284.00.30

| | |
|-------------|--------------|
| SEAT HEIGHT | 42 cm – 57cm |
|-------------|--------------|

| | |
|-------------------|---------------|
| SEAT DEPTH, WIDTH | 40 cm x 40 cm |
|-------------------|---------------|

| | |
|----------------------------|---------------|
| REQUIRED SPACE (FOOTPRINT) | 45 cm x 53 cm |
|----------------------------|---------------|

| | |
|--------|--------|
| WEIGHT | 3.6 kg |
|--------|--------|

| | |
|-----------------|--------|
| MAX BODY WEIGHT | 130 kg |
|-----------------|--------|



Australian Government
Department of Health
Therapeutic Goods Administration

ARTG 297301

info@rebotec.com.au
rebotec.com.au
218/354 Eastern Valley Way
Chatswood 2067 NSW

