

REBOTEC

GENF

24" wheel, self propelled shower commode



- Designed for easy independent manoeuvrability with well positioned brake handles, height adjustable, swinging & detachable foot rests
- Optimal over toilet, beside and shower use. Also doubles as a wheelchair
- Swinging armrests allow for ease of transfers
- Includes padded toilet seat with padded removable seat cover
- Includes removable pail and rail mechanism

MPN

Red
346.44.10
346.30.10

Blue
346.54.10

ACCESSORIES

Bucket Pail
301.00

Safety Belt
322.30.10

Armrest Cushion
440.10.32

Commode PU seat without hygiene opening
440.10.24

WIDTH, DEPTH 71cm, 95cm

SEAT WIDTH, DEPTH, HEIGHT 46cm, 45cm, 55cm

TOILET HEIGHT WITHOUT BUCKET RAILS 47cm

TOILET HEIGHT WITH BUCKET RAILS 43cm

NET WEIGHT 24kg

MAX BODY WEIGHT 130kg



Australian Government
Department of Health
Therapeutic Goods Administration

ARTG 297301

info@rebotec.com.au
rebotec.com.au
218/354 Eastern Valley Way
Chatswood 2067 NSW

