

# REBOTEC

## HAMBURG

Height adjustable, mobile commode chair



- Ergonomically shaped backrest with integrated push handle
- Easy lateral transfer through swing-away armrests
- Toilet seat with seat cushion
- Removable bucket rails
- Bucket pail included
- Height-adjustable, swing-away and removable footrests
- 5-inch swivel wheels, rear with parking brake and directional lock
- Easy to clean and disinfect

### MPN

Blue  
356.54.00

Red  
356.44.00

### ACCESSORIES

Bucket Pail  
301.00

Safety Belt  
322.30.10

IV Pole  
303.62

Armrest Cushion  
440.10.32

Commode PU seat without hygiene opening  
440.10.24

WIDTH, DEPTH 56cm, 87cm

SEAT WIDTH, DEPTH 46cm, 48cm

SEAT HEIGHT 54-59cm

TOILET HEIGHT WITHOUT BUCKET RAILS 48-53cm

NET WEIGHT 14.5kg

MAX BODY WEIGHT 130kg



Australian Government  
Department of Health  
Therapeutic Goods Administration

ARTG 297301

info@rebotec.com.au  
rebotec.com.au  
218/354 Eastern Valley Way  
Chatswood 2067 NSW

