

REBOTEC

KIEL

Height adjustable commode chair



- Designed for over toilet, in shower and beside use
- Folding armrests accommodate ease of side transfers
- Includes padded removable seat cushion
- Includes removable pail and rail device
- Height adjustable, swinging and detachable foot rests
- Sturdy rubber ferrules for safe standing

MPN

Blue
347.54.00

ACCESSORIES

Bucket Pail
301.00

Safety Belt
322.30.10

IV Pole
303.62

Armrest Cushion
440.10.32

Commode PU seat without hygiene opening
440.10.24

WIDTH, DEPTH	57cm, 63cm
--------------	------------

SEAT WIDTH, DEPTH	46cm, 45cm
-------------------	------------

SEAT HEIGHT	51-56cm
-------------	---------

NET WEIGHT	11.2kg
------------	--------

MAX BODY WEIGHT	130kg
-----------------	-------



Australian Government
Department of Health
Therapeutic Goods Administration

ARTG 297301

info@rebotec.com.au
rebotec.com.au
218/354 Eastern Valley Way
Chatswood 2067 NSW

