## REBOTEC KOLN

## Commode chair



- Swing open armrests for easy side transfer
- Includes seat cover cushion
- Includes removable pail
- Designed for use at bedside, in the shower and as a toilet support

## MPN

Blue 340.54.00

Red 340.44.00

## **ACCESSORIES**

Bucket Pail 301.00

Safety Belt 322.30.10

IV Pole 303.62

Armrest Cushion 440.10.32

Commode PU seat without hygiene opening 440.10.24

| WIDTH, DEPTH      | 57cm, 63cm |
|-------------------|------------|
| SEAT WIDTH, DEPTH | 46cm, 45cm |
| SEAT HEIGHT       | 52cm       |
| NET WEIGHT        | 11.2kg     |
| MAX BODY WEIGHT   | 130kg      |





info@rebotec.com.au rebotec.com.au 218/354 Eastern Valley Way Chatswood 2067 NSW

