

# REBOTEC

## MONACO

### Shower stool with hygiene opening



- Height adjustable and small foot print for limited space
- Easy tool free assembly for travelling
- Front cut out for easy personal hygiene access

#### MPN

Monaco  
363.00.10

#### ACCESSORIES

Chair Travel Bag  
284.00.30

SEAT HEIGHT	42 cm – 57cm
-------------	--------------

SEAT DEPTH, WIDTH	40 cm x 40 cm
-------------------	---------------

REQUIRED SPACE (FOOTPRINT)	45 cm x 45 cm
----------------------------	---------------

WEIGHT	2.8 kg
--------	--------

MAX BODY WEIGHT	130 kg
-----------------	--------



Australian Government  
Department of Health  
Therapeutic Goods Administration

ARTG 297301

info@rebotec.com.au  
rebotec.com.au  
218/354 Eastern Valley Way  
Chatswood 2067 NSW

