

REBOTEC

SAFE-IN-KIDS

Forearm cuff and hinge crutches for kids



- Forearm cuff and hinge for kids
- Ergonomic handle with reflector
- Bright, fun and colourful design for kids
- 10 height and 3 handle to cuff length adjustments

MPN

Blue
122.95.10

Yellow
122.96.10

ACCESSORIES

Upholstered Crutch Handle Pad
160.01.01

Blue 17mm Ferrules
200.50.10

Yellow 17mm Ferrules
200.20.10

Black 17mm Ferrules
200.10.10

Crutch Strap
167.00.16
167.00.30
167.00.50

HEIGHT FROM HANDLE TO CUFF 16-21cm

HEIGHT OF HANDLE 50-78cm

DIAMETER Ø IN MM 20/16.9

MAX BODY WEIGHT 100 kg



Australian Government
Department of Health
Therapeutic Goods Administration

ARTG 297302

info@rebotec.com.au
rebotec.com.au
218/354 Eastern Valley Way
Chatswood 2067 NSW

