

# REBOTEC

## VENEDIG

### Self-propelled commode for small adults



Self propelled commode chair for smaller adults and teens. with narrower commode aperture and softer backrest than regular self propelled chairs.

- PU cushioned seat and backrest, no cracking or reupholstering
- Dignity hygiene aperture cover included
- Commode pail bucket included
- Strong integrated armrests for side transfer
- Detachable foot support

#### MPN

Blue  
380.54.10

#### ACCESSORIES

Bucket Pail  
301.00

Safety Belt  
322.30.10

Armrest Cushion  
440.10.32

Commode PU seat without hygiene opening  
440.10.24

WIDTH, DEPTH 67cm, 97cm

SEAT WIDTH, DEPTH, HEIGHT 42cm, 44cm, 55cm

MAX. TOILET WIDTH 33cm

TOILET HEIGHT WITH BUCKET RAILS 43cm

NET WEIGHT 23.4kg

MAX BODY WEIGHT 100kg



Australian Government  
Department of Health  
Therapeutic Goods Administration

ARTG 297301

info@rebotec.com.au  
rebotec.com.au  
218/354 Eastern Valley Way  
Chatswood 2067 NSW

