

REBOTEC

BOSTON

Mobile shower chair, 150kg rated



- Ergonomic backrest with backrest with carer grip
- Swinging armrests for easy side transfer
- Seat with included PU cushion
- Included, removable bucket rails removable
- Bucket pail included
- Height adjustable, swinging out of the way and detachable foot support
- 4 swivelling castors 5-inch with parking/direction brake

MPN

Blue
342.54.00

Grey
342.03.00

ACCESSORIES

Bucket Pail
301.00

Safety Belt
322.30.10

IV Pole
303.62

Armrest Cushion
440.10.32

Commode PU seat without hygiene opening
440.10.24

WIDTH, DEPTH 62cm, 89cm

SEAT WIDTH, DEPTH, HEIGHT 52cm, 48cm, 53cm

TOILET HEIGHT WITHOUT BUCKET RAILS 47cm

TOILET HEIGHT WITH BUCKET RAILS 44cm

NET WEIGHT 16.2kg

MAX BODY WEIGHT 150kg



Australian Government
Department of Health
Therapeutic Goods Administration

ARTG 297301

info@rebotec.com.au
rebotec.com.au
218/354 Eastern Valley Way
Chatswood 2067 NSW

