

REBOTEC

NEW WALK

Crutches with shock absorbing springs



- Crutches with spring shock absorbers
- 10 height adjustments and 3 cuff length adjustments
- Internal shock absorbing and dampening spring
- 2 spring damper settings, no tools required

MPN

Blue Ergonomic
125.50

ACCESSORIES

Upholstered Crutch Handle Pad
160.01.02

Blue 19mm Ferrules
200.50.00

Big Foot Ferrules
168.10.17

Blue Crutch Lanyard
167.00.50

HEIGHT FROM HANDLE TO CUFF 22-25cm

HEIGHT OF HANDLE 80-100cm

DIAMETER Ø IN MM 22/19

MAX BODY WEIGHT 130 kg



Australian Government
Department of Health
Therapeutic Goods Administration

ARTG 297302

info@rebotec.com.au
rebotec.com.au
218/354 Eastern Valley Way
Chatswood 2067 NSW

