

# REBOTEC

## PHOENIX MULTI

### Electric Tilt-in-Lift power shower chair







Rebotec Phoenix Multi in standard configuration includes seat with hygiene opening, adjustable headrest, bucket and bucket device and a charging charging station with second, backup battery. The chair is certified class IP56 (protected from limited dust ingress, protected from high pressure water jets from any direction).

- Swing open armrests for easy side transfer
- Cushioned seat, leg and backrest with hygiene opening
- Tilt 90-130 degrees supine for wash and hoist transfer position
- Lift to eye level, prevent back injury, im prove dignity
- 150 kg certified capacity, multi award winning design

### MPN

Blue 352.54.10

#### ACCESSORIES

Bucket Pail 301.00

Safety Belt 322.30.10

IV Pole 303.62

Armrest Cushion 440.10.32

WIDTH, DEPTH 68cm,112cm

SEAT WIDTH, DEPTH 54cm, 43cm
SEAT HEIGHT 60-100cm

RATING IP56

NET WEIGHT 53.6kg

MAX BODY WEIGHT 150kg





info@rebotec.com.au rebotec.com.au 218/354 Eastern Valley Way Chatswood 2067 NSW

