


REBOTEC

QUICK'N EASY

Underarm crutches



- Comfortable cushioned armpit support
- Height adjustable pushpin base & clip lock hand-grips
- Soft and stable comfort rubber hand-grips
- Designed to accommodate weights of various users up to 130kg
- Sold in pairs

MPN

Child
127.10.10x2

Youth
127.20.10x2

Adult
127.30.10x2

Tall
127.40.10x2

ACCESSORIES

Black 19mm Ferrules
200.10.00

Red 19mm Ferrules
200.30.00

Aqua 19mm Ferrules
200.80.00

Big Foot Ferrules
168.10.17



Australian Government
Department of Health
Therapeutic Goods Administration

ARTG 297302

info@rebotec.com.au
rebotec.com.au
218/354 Eastern Valley Way
Chatswood 2067 NSW



Child Youth Adult Tall

SHOULDER 83-103cm 101-121cm 118-138cm 136-156cm

HANDLE 48-71cm 58-81cm 68-91cm 75-98cm

DIAMETER Ø 19/22 19/22 19/22 19/22