

# REBOTEC

## SAFE-IN-SOFT

### Forearm crutches with cuff & hinge



- Hinged Forearm cuff for improved mobility
- 10 adjustments in height + 3 height adjustments from handle to cuff
- Safety night light reflector facing back
- 130 kg weight capacity

#### MPN

Light Grey  
117.00

Black  
117.10

Red  
117.30

Blue  
117.50

Aqua  
117.80

#### ACCESSORIES

Upholstered Crutch Handle Pad  
160.01.02

19mm Ferrules  
200.00.xx

Big Foot Ferrules  
168.10.17

Crutch Lanyard  
167.00.xx

HEIGHT FROM HANDLE TO CUFF 22-25 cm

HEIGHT OF HANDLE 75-97 cm

DIAMETER Ø IN MM 22/19

MAX BODY WEIGHT 130 kg



Australian Government  
Department of Health  
Therapeutic Goods Administration

ARTG 297302

info@rebotec.com.au  
rebotec.com.au  
218/354 Eastern Valley Way  
Chatswood 2067 NSW

