

# REBOTEC

## YANO

### Arthritic, rheumatism forearm crutches



- Wide velcro fastener forearm support
- 10 height adjustment levels
- Forearm and handle inclination adjustment
- Rotating handles
- Solid high diameter ferrule
- Sold in pairs (2pc)
- 150 kg body weight capacity

#### MPN

Grey  
130.00.10x2

#### ACCESSORIES

Big Foot Ferrules  
168.10.17

HEIGHT OF FOREARM SUPPORT	92-114cm
---------------------------	----------

HEIGHT OF HANDLE	max. 10 cm
------------------	------------

DIAMETER Ø IN MM	22/19
------------------	-------

MAX BODY WEIGHT	150 kg
-----------------	--------



Australian Government  
Department of Health  
Therapeutic Goods Administration

ARTG 297302

info@rebotec.com.au  
rebotec.com.au  
218/354 Eastern Valley Way  
Chatswood 2067 NSW

