

# REBOTEC

## SAFE-IN-EXCESS LENGTH

Tall forearm crutches with cuff & hinge



- Hinged forearm cuff for improved mobility
- 10 adjustments in height + 3 height adjustments from handle to cuff
- Safety night light reflector facing back

### MPN

Black  
116.10.20

### Accessories

Handle Pads  
160.01.02

19mm Ferrules - Black  
200.10.00

Crutch Strap - Charcoal  
167.00.16

Big Foot Ferrules - 18.5-20.5mm  
168.10.17



Australian Government  
Department of Health  
Therapeutic Goods Administration

ARTG 297302

HEIGHT FROM HANDLE TO CUFF 22-25 cm

HEIGHT OF HANDLE 85-107 cm

DIAMETER Ø IN MM 22/19

info@rebotec.com.au  
rebotec.com.au  
218/354 Eastern Valley Way  
Chatswood 2067 NSW

