

REBOTEC

SWING

Bathroom grab rail



- Swing grab rails can help make your bathroom safer and more accessible
- Rails can be mounted diagonally for support near a toilet or mounted vertically or horizontally for a shower support and safety
- 4 lengths to choose from to fit any bathroom


MPN

60s
604.05.45

75s
604.05.60

90s
604.05.75

105s
604.05.90



60s Length A: 45 cm Length B: 67 cm

75s Length A: 60 cm Length B: 82 cm

90s Length A: 75 cm Length B: 97 cm

105s Length A: 90 cm Length B: 112 cm



Australian Government
Department of Health
Therapeutic Goods Administration

ARTG 305316

info@rebotec.com.au
rebotec.com.au
218/354 Eastern Valley Way
Chatswood 2067 NSW

