

REBOTEC

CUFF UPHOLSTERY

Upholstered crutch cuff pad



- Crutch cuff upholstered pad
- Pressure-reducing and skin-friendly material
- Slip-resistant neoprene with nylon cover
- Easy Velcro attachment around crutch handle
- Hand wash in warm water

MPN

Large
160.02.03x2



Australian Government
Department of Health
Therapeutic Goods Administration

ARTG 297302

info@rebotec.com.au
rebotec.com.au
218/354 Eastern Valley Way
Chatswood 2067 NSW

