

REBOTEC

SAFE-IN-ANATOM-SOFT

Crutches anatomically shaped handle



Forearm crutches with cuff & hinge

- Hinged Forearm cuff for improved mobility
10 adjustments in height + 3 height adjustments from handle to cuff
- Safety night light reflector facing back
- 130kg weight capacity

MPN

Light Grey
118.00

Black
118.10

Red
118.30

Blue
118.50

Aqua
118.80

ACCESSORIES

Upholstered Crutch Handle Pad
160.01.04

19mm Ferrules
200.00.xx

Big Foot Ferrules
168.10.17

Crutch Lanyard
167.00.16
167.00.30
167.00.50

HEIGHT FROM HANDLE TO CUFF 22-25 cm

HEIGHT OF HANDLE 75-97 cm

DIAMETER Ø IN MM 22/19

MAX BODY WEIGHT 130 kg



Australian Government
Department of Health
Therapeutic Goods Administration

ARTG 297302

info@rebotec.com.au
rebotec.com.au
218/354 Eastern Valley Way
Chatswood 2067 NSW

